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Baseline					
Collection	Manufacturer/ number	product	Handling/processing in the study Center	Pipetting procedure	Storage
<b>Blood samples</b>					
2x 10 ml serum tube with clotting activator without gel	BD Life Science, Franklin Lakes, USA / 367896		4 x inverted, clotting time 30-45 minutes, centrifugation 10 minutes, 15°C, 2500xg	Pipetting roboter (Hamilton Robotics, Reno, USA), plates with 96 wells (Azenta Life Sciences, Chelmsford, USA) 250 µl and 500 µl tubes )	<ul style="list-style-type: none"> <li>Central biorepository (2/3 of all collected samples):               <ul style="list-style-type: none"> <li>500µl at -80°C</li> <li>250µl tubes at -180°C</li> </ul> </li> <li>Local back-up storage (1/3 of all collected samples)               <ul style="list-style-type: none"> <li>-80°C/-180°C</li> </ul> </li> </ul>
3x 10 ml K <sub>2</sub> EDTA tube (1,8mg/ml)	BD Life Science, Franklin Lakes, USA / 367525		4 x inverted, 30 seconds – max. 5 minutes blood mixer, interims storage until shipping by room temperature	Pipetting roboter, plates with 96 wells, separated into: <ul style="list-style-type: none"> <li>Plasma (250 µl)</li> <li>Erythrocytes (250 µl tubes)</li> <li>Extended buffy coat (for later DNA isolation) (500 µl tubes)</li> </ul> <hr/> <ul style="list-style-type: none"> <li>DNA isolated from 1 buffy coat (expected to be completed in September 2022)</li> </ul>	<ul style="list-style-type: none"> <li>Central biorepository (2/3 of all collected samples)               <ul style="list-style-type: none"> <li>-180°C</li> </ul> </li> <li>Local back-up storage (1/3 of all collected samples)               <ul style="list-style-type: none"> <li>-80°C/-180°C</li> </ul> </li> <li>Central biorepository (2/3 of all collected samples)               <ul style="list-style-type: none"> <li>-180°C</li> </ul> </li> <li>Local back-up storage (1/3 of all collected samples)               <ul style="list-style-type: none"> <li>-80°C/-180°C</li> </ul> </li> <li>Central biorepository (all collected samples)               <ul style="list-style-type: none"> <li>2x500µl -80°C, one used for DNA isolation</li> <li>4x500µl -180°C</li> <li></li> </ul> </li> <li>Central biorepository               <ul style="list-style-type: none"> <li>1 x 250/400 µl</li> </ul> </li> </ul>
2.5 ml RNA stabilized blood with PAXgene	BD Life Science, Franklin Lakes, USA / 762165		10 x inverted, 1-2 hours stored at room temperature, 24 hours -max. 72 hours by -20°C, then put into -80°C		<ul style="list-style-type: none"> <li>central biorepository (all collected samples) -80°C (semi-automated storage system) Prior to storage, the RNA samples were volume reduced in central biorepository. The workflow is as follows: thaw, transfer to 50 ml flacons, centrifuge, remove supernatant, collect in resuspension buffer and vortex, transfer to FluidX tubes</li> </ul>
<b>Urine samples</b>					
100 ml collecting cup	BD Life Science, Franklin Lakes, USA / Urine collection cup 364941, Urine tube / 365000		10 ml transferred into urine vacutainer, centrifugation 10 minutes, 15 °C, 2500xg	Pipetting roboter, plates with 96 wells (250 µl and 500 µl tubes)	<ul style="list-style-type: none"> <li>central biorepository (2/3 of all collected samples)               <ul style="list-style-type: none"> <li>500µl at -180°C</li> <li>250µl at -180°C</li> </ul> </li> <li>local back-up storage (1/3 of all collected samples)               <ul style="list-style-type: none"> <li>-80°C/-180°C</li> </ul> </li> </ul>

Saliva samples			
Paraffin chewing gum	GC Germany, Bad Homburg, Germany,	Participant chews gum for one minute, then saliva is dispensed; if necessary, the chewing time is extended by another minute	Pipetted by hand, plates with 96 wells (500 µl tubes) <ul style="list-style-type: none"> <li>• central biorepository (all collected samples) <ul style="list-style-type: none"> <li>○ -80°C</li> </ul> </li> </ul>
Medicine cup with lid	Sarstedt, Nümbrecht, Germany / 75.1337.500 + 76.1340.560		
Nasal swabs			
Sterile nylon flock swabs	Copan Flock Technologies, Brescia, Italy	Two swabs are transferred in one 2 ml cryotube pre-filled with RNA-later	- <ul style="list-style-type: none"> <li>• central biorepository (all collected samples) <ul style="list-style-type: none"> <li>○ -80°C</li> </ul> </li> </ul>
Stool samples			
At the same day or one day before the visit in the study center, at home with stool catcher for sample collection	Süsse Labortechnik, Gudensberg, Germany / S1000	Stabilized stool sample is sent to the central university laboratory Greifswald, where it is aliquoted in 500 µl tubes	- <ul style="list-style-type: none"> <li>• central biorepository (all collected samples) <ul style="list-style-type: none"> <li>○ -80°C</li> </ul> </li> </ul>
Two stool collecting tubes were filled: <ul style="list-style-type: none"> <li>• one is without additives</li> </ul>	Biosepar, Mühldorf am Inn, Germany		
<ul style="list-style-type: none"> <li>• one tube was prefilled with 3.5ml RNA-later (RNA stabilizing additive with K2EDTA)</li> </ul>	Quiagen, Hilden, Germany		
The samples are packed in a protective bag "Süsse Schutzbeutel" in a box "Süsse Post Box Slim" and brought to the study center at the time of the examination	Süsse Labortechnik, Gudensberg, Germany / H64490 and H11227		

Second examination					
Collection	Manufacturer/ number	product	Handling/processing in the study Center	Pipetting procedure	Storage
<b>Blood samples</b>					
1x 10 ml serum tube with clotting activator without gel	BD Life Science, Franklin Lakes, USA / 367896		4 x inverted, clotting time 30-45 minutes, centrifugation 10 minutes, 15°C, 2500xg	Pipetting roboter (Hamilton Robotics, Reno, USA), plates with 96 wells (Azena Life Sciences, Chelmsford, USA) 250 µl tubes)	<ul style="list-style-type: none"> <li>• Central biorepository (2/3 of all collected samples):               <ul style="list-style-type: none"> <li>○ 1x250µl tube at -80°C*</li> <li>○ 7x250µl tubes at -180°C</li> </ul> </li> <li>• Local back-up storage (1/3 of all collected samples)               <ul style="list-style-type: none"> <li>○ -80°C/-180°C</li> </ul> </li> </ul>
2x 10 ml K <sub>2</sub> EDTA tube (1,8mg/ml)	BD Life Science, Franklin Lakes, USA / 367525		4 x inverted, 30 seconds – max. 5 minutes blood mixer, interims storage until shipping by room temperature	Pipetting roboter, plates with 96 wells, separated into: <ul style="list-style-type: none"> <li>• Plasma (250 µl tubes)</li> <li>• Extended buffy coat (for later DNA isolation) (500 µl tubes)</li> </ul>	<ul style="list-style-type: none"> <li>• Central biorepository (2/3 of all collected samples)               <ul style="list-style-type: none"> <li>○ 250 µl at -180°C</li> </ul> </li> <li>• Local back-up storage (1/3 of all collected samples)               <ul style="list-style-type: none"> <li>○ -80°C/-180°C</li> </ul> </li> <li>• Central biorepository (2/3 of all collected samples)               <ul style="list-style-type: none"> <li>○ 1x500µl -80°C**</li> <li>○ 1x500µl -180°C</li> </ul> </li> <li>• Local back-up storage (1/3 of all collected samples)               <ul style="list-style-type: none"> <li>○ -80°C/-180°C</li> </ul> </li> </ul>
<b>Urine samples</b>					
100 ml collecting cup	BD Life Science, Franklin Lakes, USA / Urine collection cup 364941, Urine tube / 365000		10 ml transferred into urine vacutainer, centrifugation 10 minutes, 15 °C, 2500xg	Pipetting roboter, plates with 96 wells (250 µl)	<ul style="list-style-type: none"> <li>• central biorepository (2/3 of all collected samples)               <ul style="list-style-type: none"> <li>○ 250µl at -180°C</li> </ul> </li> <li>• local back-up storage (1/3 of all collected samples)               <ul style="list-style-type: none"> <li>○ -80°C/-180°C</li> </ul> </li> </ul>

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Sample type	Volume	Aliquots per subject (Total/in Biorep)	Study participants (Completeness, %)	Study participants at central storage (Completeness, %)	Retrievable	Preprocessing required	Preferred analysis
Baseline examination*							
Plasma	250 µl	48/32	95.7	94.5	Yes (-180°C storage temperature)	No	
Serum	250 µl	30/20	96.8	95.6	Yes (-180°C storage temperature)	No	
Serum	500 µl	1/1	95.7	94.1	Yes (-180°C storage temperature)	No	
Erythrocytes	250 µl	6/4	95.5	94.3	Yes (-180°C storage temperature)	No	
Buffy Coat	500 µl	6/6	94.9	93.5	Yes (-80°C and -180°C storage temperature)	DNA isolation, extraction of DNA is currently underway at the HMGU	DNA isolation → genomics
Urine	250 µl	12/8	96.7	95.5	Yes (-180°C storage temperature)	No	
Urine	500 µl	4/4	95.9	94.4	Yes (-180°C storage temperature)	No	
Saliva	500 µl	2/2	93.5	92.2	Yes (-80°C storage temperature)	DNA isolation, to be done by user	DNA isolation → genomics (microbiome)
Nasal Swab	2 ml	1/1	95,7	93,6	Yes (-80°C storage temperature)	DNA isolation, to be done by user	DNA isolation → genomics (microbiome)
Native Stool	5 ml	1/1	75.7	75.3	Yes (-80°C storage temperature)	No	
Stabilised Stool	500 µl	4/4	72.5	71.9	Yes (-80°C storage temperature)	RNA isolation	RNA isolation → transcriptomics
RNA (PAX Gene)	10 ml	1/1	21.3	19.1	No, will be volume reduced (500 µl) at the central Biorepository	-	-
RNA (500 µl)	500 µl	1/1	72.1	72.1	Yes, depending on the state of the RNA volume reduction (-80°C storage temperature)	RNA isolation, to be done by user	RNA isolation → transcriptomics
DNA	500 µl	1/1	41,2	41,2	Depending on the state of DNA isolation from Buffy Coat samples		Any genomic analysis

Follow-up examination**							
Plasma	250 µl	18/12	83.3	71,6	Yes (-180°C storage temperature)	No	
Serum	250 µl	12/8	83.3	71.6	Yes (-80°C samples since August 2022 and -180°C storage temperature)	No	
Buffy Coat	500 µl	2/2	82.9	69.7	Yes (-180°C storage temperature)	DNA isolation, extraction will take place at HMGU, stored at -80 C	DNA isolation → genomics
Urine	250 µl	6/4	84.7	72.9	Yes (-180°C storage temperature)	No	

Numbers are subject to change when samples are released from stock or consent is withdrawn.

\*Planned number of participants for baseline examination has been set to 205,000 except for stool samples. Collection of stool samples had not been introduced by the start of the study, but only in October 2016. In the case of the stool samples, the planned number of participants was set at 100,000.

\*\*Planned number of participants for follow-up examination has been set to 135,000.

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